\*FattaLab<sup>®</sup>available for sale in Hong Kong only.

## LIGHTWEIGHT PORTABLE PERSONALIZED



# FATTY LIVER DIAGNOSTIC DEVICE

MODEL: FL-CC M1

#### THE WORLD'S FIRST LIGHTWEIGHT FATTY LIVER INTELLIGENT ASSESSMENT SYSTEM

PIONEERING THE ERA OF HEALTH MONITORING 3.0

Smart Weight Tracking 1.0 Precision Fat Mapping 2.0 Visceral Fat Quantification 3.0

Lab precision,

Home convenience.

# **CONTACT US**

#### **EIELING TECHNOLOGY LIMITED (Headquarters)**

Room 529, 5/F, Core Building 5W, Hong Kong Science Park, New Territories, Hong Kong +852 5562 1051 info@eieling.com www.eieling.com

Hong Kong | Shanghai | Shenzhen | Nanjing

 $\bigtriangledown$ 

Get quantitative liver fat results in 30 seconds

One-button measurement with automated analysis and reports via Smart App

### CLINICAL APPLICATIONS

#### Medical Centers | Clinics

FattaLab® Integrate Testina into Community Screening Programs. Targeted at high-risk populations (diabetes, obesity, etc.) for fatty liver disease. program correlates the FattaLab® data with key metabolic markers (blood glucose, lipids, etc.) to build a predictive metabolic syndrome model. This enables early identification of at-risk patients, while centers provide continuous intervention strategies through monthly monitoring reports.

#### **Fitness Centers**

Launch a Fatty Liver Reversal Program and offer members a "Liver Health Rewards System" to improve retention rates.

#### **Insurance Providers**

Integrate FattaLab<sup>®</sup> testing into family health insurance plans to lower user costs, Insurance providers may use premium pricing models that reduce fees for low-risk clients depending on FattaLab<sup>®</sup>-generated data.

Empowering everyone with data-driven precision health management.



From 'Estimation' to 'Diagnostic-grade' Precision

Auto Intra-Focus Fusion (AIFF) Human Liver Fat Detection Algorithm

\*The model was trained on real-world datasets encompassing over 10,000 clinical cases.



# **KEY ADVANTAGES**



The world's first technology to combine time-domain and frequency-domain modulation, reducing subcutaneous fat interference for obese or high body fat users.



Al-powered liver region recognition focuses on hepatic parenchyma while excluding non-hepatic regions.



Multi-step fine-tuning mechanism overcomes small right liver lobe limitations in large-ROI acquisitions.

### PERFORMANCE PARAMETERS

\*Refer to the product manual for definitive specifications.

TECHNICAL SPECIFICATIONS	PARAMETERS
Dimensions	130*58*25 mm
Net Weight	140 g
Maximum Detection Depth	200 mm
Measurement Depth	20~100 mm
Measurement Range	100-400 dB/m
Battery life	≥3h